mgMed basic

Terapeutico

CE

1370

ITA - ENG - FRA - ESP - DEU

HSD

ITALIANO	Manuale d'uso e programmi per l'uso terapeutico	PAGE 2-16
ENGLISH	User manual and programs for therapeutic use	PAGE 17-30
FRANÇAIS	Manuel d'utilisation et programmes à usage thérapeutique	PAGE 31-44
ESPAÑOL	Manual de usuario y programas para uso terapéutico	PAGE 45-58
DEUTSCH	Gebrauchsanweisung und Programme für den therapeutischen Einsatz	PAGE 59-72

ITALIANO

Fabbricante:

HSD S.r.l.

Strada Acquasalata, 11 - int. 6 - 47899 Serravalle RSM

Tel (+378) 0549.941496 - Fax (+378) 0549.976456

Numero Verde : 800 02 00 55

Mail : commerciale@hds.sm

Web : www.hsd.sm

Gentile Cliente,

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Dear Customer,

we congratulate you on your choice.

It now has at its disposal a high-performance product, built according to the most recent and severe European standards, and capable of guaranteeing a high standard of quality and safety.

Please pay close attention to this manual by reading it in its entirety in order to familiarize yourself with the appliance.

HSD S.r.l.

Distributed by:

Introduction

This technical document describes the rules for the correct use and user maintenance of the appliance.

Compliance with the rules and instructions given here will allow it to be used in the safest possible way.

Information

If, after reading, doubts or uncertainties persist about the operations to be carried out, it is necessary to contact the trusted Reseller or the Manufacturer.

No part of this manual may be reproduced, copied or disclosed by any means without the written authorization of HSD S.r.l. Any infringements will be prosecuted in the manner and within the times provided for by the law. All names and trademarks mentioned in this manual are the property of their respective manufacturers.

General index

Description	Page
Introduction	17
General index	17
Description of the components	18
Warranty	19
Symbolology used	19
General safety information	19
Contraindications	19
Technical and functional characteristics	19
Technical data	20
Identification	20
Put in action	20
Language setting procedure	21
Maintenance	22
Accessories and spare parts	23
Ultrasound therapy	23
What is cellulite?	24
Programs	25
Questions and answers	29
Troubleshooting	30
Reference regulations	30
Declaration of conformity	30

Description of the components

- | | | | |
|---|---|----|---|
| 1 | DISPLAY - used for dialogue with the operator. | 7 | Led "ACTIVE" - indicates the appliance is in |
| 2 | "START / STOP" button - used to start, pause and stop the treatment. | 8 | operation. Treatment head connector |
| 3 | "PROGRAMS + , PROGRAMS -" buttons - allow you to move between programs | 9 | Treatment head |
| 4 | Led "ON" - indicates that the appliance is connected to the electrical network but not in operation. | 10 | EDIT button (useful for changing language) |
| 5 | Led "STOP" - indicates that the appliance is not in operation. | 11 | Power connector housing |
| 6 | Led "PAUSE" - indicates that the appliance is in operation but temporarily stopped. | 12 | General ignition switch. |

Package Contents

- Appliance 220 ÷ 240 VAC 50 - 60 Hz
- Optional accessories based on configuration
- Power cable
- User manual with treatment guide and warranty certificate
- Transport bag

The device is supplied with all the equipment indicated above, so once the packaging has been opened, check that the equipment is complete. If any element is missing, contact the Seller.

OPTIONAL

50mm head

OPTIONAL

40mm head

OPTIONAL

40mm tip head

OPTIONAL

Double head diffuser

⚠ Caution!

Dual head diffusers should never be used in immersion programs

Warranty

Each device is subject to careful checks and is tested before being put on sale. The consumer is the owner of the rights specified by current legislation. The warranty has a duration of two years from the date of purchase and is extended to all parts, except those subject to wear. It lapses if the HSD device is not used appropriately or tampered with by unauthorized personnel, or if the defect is not reported within two months from the date on which it appeared. The parts deemed by HSD S.r.l. defective at the time of manufacture will be replaced free of charge.

HSD S.r.l.

Useful Life

HSD S.r.l. provides that the duration of maintenance of the functionality and safety of the device in the conditions of use, maintenance care indicated in the User Manual is 6 years.

The devices are guaranteed, with correct use and maintenance, by the HSD S.r.l. for 24 months.

HSD S.r.l.

Symbology used

Caution

Draws attention to situations or problems that could affect safety.

Information

Draws attention to important situations of a general nature.

General safety information

Some general safety guidelines are summarized below:

- carefully read the instructions for use before starting the appliance;
- do not use the appliance in the presence of flammable mixtures with air or with oxygen or with nitrous oxide, absolutely avoid explosive environments and humid environments or environments in which large quantities of water are present, as the apparatus is not protected against the abundant penetration of liquids.
- check, before each treatment, the integrity of the power cable, the treatment head cable and connectors;
- to power the device, use only the supplied bipolar cable or request it from HSD S.r.l.;
- never carry out any operation on the appliance with damp or wet hands;
- do not immerse the appliance in water for cleaning;
- consult the paragraph on the available programs before using the device;
- do not place weights, and in particular containers with liquid on the appliance, do not spill liquids on it;
- in case of abnormal behavior or malfunction, stop the treatment and turn off the device, do not try to open or repair the device. Contact HSD S.r.l. for any type of intervention or control.

WEEE: This product complies with Directive 2012/19 / EU.

The crossed-out bin symbol shown on the appliance indicates that the product, at the end of its useful life, having to be treated separately from household waste, must be disposed of in a separate collection center for electrical and electronic equipment or returned to the dealer when purchasing a new equivalent appliance. The user is responsible for delivering the appliance at the end of its life to the appropriate collection facilities, under penalty of penalties provided for by current waste legislation. Adequate separate collection for the subsequent forwarding of the decommissioned appliance for recycling, treatment and environmentally compatible disposal helps to avoid possible negative effects on the environment and health and favors the recycling of the materials referred to the product is composed. For more detailed information regarding the collection systems available, contact the local waste disposal service. Manufacturers and importers comply with their responsibility for environmentally compatible recycling, treatment and disposal both directly and by participating in a collective system.

Contraindications

Caution

The appliance is locally contraindicated:

- in the pregnant woman;
- in those with cardiac pacemakers;
- on tumor areas;

The appliance is not recommended:

- in the severely debilitated patient who is unable to communicate;

Caution

In case of doubts or concerns it is necessary to consult a doctor.

Intended use and applications

Ultrasound therapy is a physiotherapy technique that uses the biological effects produced by ultrasound for both therapeutic and aesthetic purposes.

Ultrasound: Frequency 1MHz, 3MHz / Power: 0.25 - 2.50W / cm² / DutyCycle: 10% -100%.

The device is studied and designed to be used in environments such as:

1. Outpatient 2. Physiotherapeutic 3. Aesthetic 4. Domestic

The use of this device is allowed to medical, paramedical and to the patient himself if suitably informed of the methodologies used by the unit. In case of doubts or concerns, it is necessary to consult a doctor.

Specific therapeutic applications

- Neck pain of arthritic origin or tension muscle; - Cervical brachial neuralgia, lumbosacral pains from arthrosis or from minor intervertebral disorders;
- Low back pain of non-discal origin; - arthrosis of the hips; - arthrosis of the knees; - meniscopathies; - joint pains of various kinds; - rheumatoid arthritis of the wrists and hands; - ankylopoietic spondylitis; - scapula humeral periarthritis;
- epicondylitis-epitrocleitis; - distortions; - muscle distractions; - bruises;
- hematomas; - reflex algodystrophies; - tenovaginitis; - bunions; - snap fingers;
- painful sequelae of already welded fracture; - consolidation delays; - achillodine;
- ischial pain; - algic dorsal apophyses; - coccygodynias; - groin pain;
- talalgia; - metatarsalgia.

Specific fitness and beauty applications

- Absorption of active substances on the face and body;
- treatment of superficial wrinkles and sensitive skin;
- acne treatment; - boils on the face and body;
- anti-cellulite treatments

Caution

The device must not be used for:

1. electrodiagnosis
2. experimental use;
3. operations not included in the examples of use given in the instructions for use and maintenance.

Any use other than those described or the use of products that do not comply with HSD specifications is to be considered improper.

Technical and functional characteristics

Fabbricante: HSD S.r.l.
Strada Acquasalata, 11 - int. 6 - 47899 Serravalle RSM
T. (+378) 0549.941496 - F. (+378) 0549.976456

Technical data

POWER SUPPLY	220 ÷ 240 VAC
FREQUENCY POWER SUPPLY	50 - 60 Hz
ABSORBED POWER	50 W
ELECTRIC INSULATION CLASS	II
TYPE OF PARTS APPLIED	BF
WORKING FREQUENCY ULTRASOUND	1 - 3 Mhz
MAXIMUM OUTPUT POWER	2,5 W/cm ²
POWER REGULATION	0 ÷ 100%
MODULATION FREQUENCY	10 ÷ 100%
MODULATION TYPE	Pulsed, continue
DIAMETER OF THE EMITTER SURFACE	50 mm / 40 mm
ACTUAL RADIATION AREA TIMER MAX	6,15 cm ² / 3,79 cm ²
THERAPEUTIC PROGRAMS	12 minuti
	51

DIMENSIONS AND WEIGHT

HEIGHT	70	mm
DEPTH	210	mm
WIDTH	320	mm
WEIGHT	3	kg

ENVIRONMENTAL CONDITIONS ALLOWED FOR OPERATION

TEMPERATURE [T]	+10 ÷ +40 °C
HUMIDITY [Ur]	30 ÷ 75 %
ATMOSPHERIC PRESSURE [mbar]	da 700 a 1060 mbar

ENVIRONMENTAL CONDITIONS ALLOWED FOR TRANSPORT AND STORAGE

For the transport and storage of the appliance the following environmental conditions are recommended:


TEMPERATURE [T]	0 ÷ 40 °C
HUMIDITY [Ur]	25 ÷ 75 %

On the back of the device and on the ultrasound head, identification plates are affixed with the technical data to be communicated to the Authorized Service Center, the Manufacturer or the Dealer in the event of a request for information or technical assistance.

Identification

- A Maker.
- B CE marking.
- C Attention, consult the attached documentation.
- D Physiological effects: elevated temperature.
- E Consult the Instructions for Use.
- F Device class II.
- G BF Type of applied parts
- H Technical data.
- I Model, serial number and year.
- L Waste disposal: treat separately from household waste.

DEVICE

HEAD 50 MM

HEAD 40 MM

DOUBLE HEAD DIFFUSER

HEAD 40 MM TIP

Information

Damage, removal or alteration of the identification plate must be immediately reported to the Manufacturer, the Authorized Service Center or the Dealer. Warranty interventions will not be possible if the identification plate is damaged, removed or altered.

Put in action

Remove the appliance from the packaging and place it on a stable horizontal rigid surface, suitable for supporting its weight, near a power outlet.

The space around the appliance must not be less than 20 centimeters on all sides.

Caution

The installation of the device is prohibited:

- in humid environments;
- near large quantities of water;
- near heat sources;
- in direct sunlight.

CONNECTION OF THE TREATMENT HEAD

Insert the connector of the treatment head into the socket provided on the front panel paying attention to the insertion guide.

It is forbidden to use the device in atmospheres that may contain flammable and / or anesthetic gases.

Caution

Never connect and disconnect the treatment head when the device is on.

CONNECTION TO THE ELECTRICITY NETWORK

Caution

Do not use multiple plugs or power adapters.

Connect the appliance to a tested and efficient socket, equipped with earth.


If in doubt, do not connect the device.

Check that the mains voltage corresponds to the indicated voltage and frequency

on the nameplate of the appliance (an incorrect voltage of the power supply can damage it without consequences for the user). The use of the appliance not properly connected to the mains it could be dangerous.

Before inserting the power supply into the socket current make sure that the n-ignition switch is positioned on "0" (zero).

Insert the power cable in housing "11" near the switch "12", making sure that it is fully inserted. Insert the connection plug into the socket.

Caution

To avoid strangulation and suffocation with the power cord: mount the unit and the power cord away from the face and out of the reach of children.

SWITCHING ON THE INSTRUMENT

Turn on the instrument by positioning the power switch on "1", the ON LED will light up to certify the correct electrical connection. Subsequently, the presentation text will scroll on the display:

						H	S	D						
			U	L	T	R	A	S	O	U	N	D	S	

and then, in the second line, the version and issue date of the operating system: E.g. : "1.00 30 Apr 09"

If the user had not connected the treatment head or had inserted it incorrectly, the display will show writings that highlight this in addition to the continuous emission of "beeps" to call attention and request the correct connection.

PREPARATION OF THE APPLIANCE

Once all the checks have been completed, a request to choose a program will be displayed

				C	H	O	O	S	E					
		A		P	R	O	G	R	A	M	M	E		

At this point, using the "▲" and "▼" keys, one of the 53 available programs will be selected and the number of the beauty program and the application point will be shown on the display:

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During the execution of the treatment program, **ACTIVE** phase, the **START / STOP** keys and the "+" and "-" keys remain enabled, with which it is possible to increase / decrease the power applied by 20% compared to the value pre-set-up.

Once the treatment program has been chosen, press **START** to start. The display will show the selected program on the first line and the treatment time, power (in%) and DutyCycle of the selected treatment on the second line. Example:

N. Progr.	Type
53	THERAPEUTIC
	TEMPORAL - MANDIBULAR PAIN

During the execution of the treatment program, **ACTIVE** phase, only the **START / STOP** button remains enabled.

After the programmed time (our 8 minutes), the machine will emit three consecutive beeps and the display will show the end of treatment message and an invitation to choose a program.

				C	H	O	O	S	E					
		A		P	R	O	G	R	A	M	M	E		

At this point, only the ON led will remain lit.

PAUSE

Pressing the **START / STOP** key switches from the **ACTIVE** phase to the **PAUSE** phase, in which the treatment is "suspended". The time will not decrease and the head will not emit signals. In this phase the **PAUSE LED** flashes while the **ACTIVE** and **START / STOP** LEDs are off. A further press of the **START / STOP** key will take us back to the **ACTIVE** phase.

STOP

Both during the **ACTIVE** phase and the **PAUSE** phase, by pressing the **START / STOP** button and holding it down for at least 1.5 seconds, I **FINISH**

definitely the treatment program.

The START / STOP led will turn off, the treatment head will not emit any signals, the duration time will rapidly decrease until it reaches zero, the "beep" will emit 3 nearby notes, the STOP led will light up and the display will show the message timed for about 3 seconds.

			P	R	O	G	R	A	M	M	E				
			E	N	D	E	D								

And below the invitation to choose a program.

				C	H	O	O	S	E				
		A		P	R	O	G	R	A	M	M	E	

At this point, only the ON led will remain lit.

Handle the treatment head with care. Careless handling can affect its characteristics.

Do not knot or twist the connection cable.

Do not use in the vicinity of the device and the treatment head, sprays, liquids or solvents.

LANGUAGE SETTING PROCEDURE

Switch on the instrument and wait until the message appears on the display:

					C	H	O	O	S	E				
		A		P	R	O	G	R	A	M	M	F		

Press and hold the **EDIT** for about 2 sec, until the display shows:

C	H	O	O	S	E		L	A	N	G	U	A	G	E	
E	N	G	L	I	S	H									

Use the "+" and "-" keys to choose between the available languages.
Press the **EDIT** button until "Choose a Program" appears in the selected language.

⚠ Caution

Handle the treatment head with care. Careless handling can affect its characteristics.
Do not knot or twist the connection cable.
Do not use in the vicinity of the device and the treatment head, sprays, liquids or solvents.

APPLICATION METHOD

The application with the appliance is simple, but in order to obtain satisfactory results it is advisable to perform the sessions taking into account some fundamental rules. We remind you that the air attenuates the ultrasounds, so it is always necessary to interpose a certain quantity of ultrasound gel or cream suitable for massage, between the treatment head and the surface of the skin to be treated. The treatment head must be applied perpendicular to the area to be treated and slid slowly, without too much pressure, collecting the gel or cream used from time to time always keeping a certain layer of gel or cream between the treatment head and the skin. To avoid overheating the treatment area it is necessary to move the head with a slow and uniform movement, taking care never to leave it stationary.

FUNDAMENTAL RULES TO RESPECT

⚠ Caution

- If it is necessary to make an "unexpected" interruption, remember to "ALWAYS" pause the machine, otherwise the head LEFT RUNNING, immobile and without creams, will overheat with the risk of damage;
- before proceeding with the application it is advisable to identify the part of the body on which to act. For anti-cellulite treatments we recommend the upright position for a correct evaluation of the area to be treated, because in a sitting or lying position there is a staggering of the silhouette;
- having identified the exact area, select the type of program and the action on the device;
- the actions that the program proposes foresee application times on small areas, therefore to act on larger areas it will take a longer amount of time, having to repeat the same program several times depending on abstention of the area to be treated;
- before starting the program, pour a large quantity of gel or a specific cream on the treatment head;
- start the program and place the treatment head on the bare part to be treated;
- move the treatment head in a slow circular motion, back and forth or as for drawing an eight. Depending on the program, you will feel a slight vibration or a slight heat.

END OF TREATMENT

At the end of the treatment:

- turn off the appliance by turning the main switch to "0" (zero). If the appliance remains on, its operation is not compromised;
- disconnect the treatment head from the device panel, taking the connector and pulling upwards perpendicularly, avoiding tearing;
- disconnect the connector plug from the socket;
- disconnect the power cable from the back of the device by pulling it outwards with very slight vertical oscillations;
- clean the treatment head, avoiding the use of solvents or abrasives, as indicated in the "Maintenance" chapter;
- collect the treatment head and power supply cables properly neat, avoiding knotting or twisting them;
- reposition the device and its accessories in the supplied case;
- store in a dry and safe place, out of the reach of children, avoiding accidental falls.

⚠ Caution

At the end of each treatment, carry out the maintenance operations described in the "Maintenance" chapter indicated with periodicity "A".

Maintenance

i Information

For maintenance, only original products and spare parts must be used. Improper use of non-original products and spare parts causes the warranty to be void.

The information in this chapter is provided in order to allow the user to obtain the best performance from the appliance in the best safety conditions.

⚠ Caution

During maintenance of the appliance:

- *never carry out any operation with damp or wet hands;*
- *do not immerse the appliance in water for cleaning, use a cloth soaked in non-abrasive neutral detergent;*
- *do not try to open or repair the appliance. Contact the Authorized Service Center for any type of intervention or check.*

MAINTENANCE OPERATION	PERIODICITY	
Check the integrity of the power cable.	M	F
Check the integrity of the treatment head cable.	M	F
Verification of the integrity of the treatment head.	A	A
Cleaning the appliance.	A.N.	A.N.
Cleaning and disinfection of the treatment head.	A	A
Periodic calibration	2Y	2Y
Electrical checks of the appliance.	2Y	2Y
REGULAR USE		
FREQUENT OR INTENSE USE		

i Information

Periodicity:

- A At each treatment.
- M Monthly.
- F Fortnightly.
- A.N. As needed.
- 2Y Biennial.

⚠ Caution

The biennial maintenance of the appliance must be carried out exclusively by HSD S.r.l.

CLEANING AND DISINFECTION

Caution

Before each action, set the ignition switch to "0" (zero) and disconnect the appliance from the power outlet.
Handle the treatment head with care. Careless handling can affect its characteristics.

CLEANING THE APPLIANCE

Remove the dust on the outside of the appliance using a soft cloth moistened with water, or with an aqueous solution consisting of 75% water and 25% denatured ethyl alcohol.

Caution

It is forbidden to use abrasives or chemicals for cleaning.

CLEANING THE TREATMENT HEAD

Disconnect the treatment head from the device and remove any residual cream or gel using a soft cloth moistened with water, or with an aqueous solution consisting of 75% water and 25% denatured ethyl alcohol. After cleaning, dry well and remove any residual moisture from the treatment head before storage or subsequent use.

Caution

It is forbidden to use abrasives or chemicals for cleaning. When cleaning, avoid contact between the connector and the cleaning solution used.

DISINFECTION OF THE TREATMENT HEAD

To disinfect the treatment head, use a cotton ball or gauze, lightly soaked in a non-irritating, alcoholic antiseptic solution (such as the type used for minor abrasions or grazes, available in pharmacies).

Caution

It is forbidden to use abrasives or chemicals for disinfection. During disinfection, avoid contact between the connector and the antiseptic solution used. Follow the directions and precautions for use provided by the manufacturer of the antiseptic solution.

CONTROL OF THE TREATMENT HEAD

At each cleaning and / or disinfection of the treatment head and in any case at the end of each treatment it is necessary to check that it does not have cracks or fissures that could allow the entry of conductive liquid. If the treatment head is cracked or cracked, replace it immediately.

Accessories and spare parts

For any orders for accessories or spare parts, consult the accessories catalog, requesting it by e-mail or by consulting the manufacturer's website.

Caution

On the back of the device and on the ultrasound head, identification plates are affixed with the technical data to be communicated to the Authorized Service Center, the Manufacturer or the Dealer in the event of a request for accessories or spare parts.

Caution

The replacement ultrasound head can only be obtained from **HSD S.r.l.**. The use of an ultrasound head other than that indicated can be dangerous.

Ultrasound therapy

Ultrasound is defined as an acoustic vibration with frequencies above the audible limit (i.e. greater than 20,000 Hz).

Ultrasonics for therapeutic purposes were introduced and systematically studied in the period following the Second World War. In Italy it spread immediately and many of our scholars have contributed significantly to experimental and clinical research in this field.

The ultrasonics spread in the form of compression-decompression waves with back and forth motion of the particles of the transmission medium, parallel to the direction of the propagation waves.

Ultrasonics are artificially produced through the piezoelectric effect using either a quartz or a disc of ceramic material. By applying electrical loads on the faces of a quartz sheet, the crystalline is compressed, reversing its direction and expanding. By subjecting the quartz to an alternating electric field, therefore, an alternation of compressions and expansions of the crystal is obtained with the production of a series of vibrations used in therapy.

When ultrasound waves travel through the tissue, they lose a certain portion of their energy: this is a process known as attenuation. The attenuation in the fabric is produced by many mechanisms: absorption, the divergence of the beam and the deflection. Absorption is the main cause of ultrasound attenuation. The ultrasound energy is absorbed by the tissue and is eventually converted into heat. For most tissues, attenuation increases with increasing frequency, so a 1.0 MHz signal would penetrate deeper than a 3.0 MHz signal, due to the lower attenuation in the tissue. Beam divergence is the degree to which the beam disperses from the transducer. The divergence of the beam decreases with increasing frequency and therefore a higher frequency signal has a more focused beam. Deflection includes the process of reflection, refraction and dispersion.

An ultrasound therapy device mainly consists of an alternating current generator (typically 1 MHz and / or 3 MHz) which supplies, via cable, an emitting treatment head in which a transducer (piezoelectric disc or aluminum foil is inserted) quartz) which converts electrical energy into mechanical energy (acoustic vibrations) which are transmitted to the tissues.

Mechanism of action

The interaction of ultrasound with biological tissues produces mechanical, thermal, chemical and cavitation effects.

Mechanical effects: the mechanical action is due to the movement of the particles of the tissues crossed by the ultrasonic wave. The pressure variations that are produced are able to determine a movement of liquids in the presence of inhomogeneity (microcurrents), an increase in membrane permeability and the disruption of the tissues by separation of the collagen fibers.

Thermal effects; the thermal effect essentially depends on two factors: the absorption characteristics of the biological medium and the reflection of energy at the interface between tissues with different acoustic impedance.

The passage of ultrasonics through the "soft" tissues creates an increase in temperature due to: - absorption linked to viscosity, - absorption due to thermal conductivity and - chemical absorption.

Chemical effects: the chemical action with modification of the local pH and permeability of cell membranes and with molecular changes is caused by the considerable acceleration forces to which the tissue particles are subjected to the passage of the ultrasonic wave.

Cavitation effects: cavitation is the ability of ultrasound to generate small dissolved gas bubbles in a fluid with subsequent increase in size and possible explosion of the bubbles. From a histological point of view, the result is an irregular cell destruction with petechial type haemorrhage. At therapeutic dosages, destructive reactions such as hemolysis would occur only in the presence of low cellular concentration and low viscosity of the medium, such as at the level of the eye and uterus.

General methods of application

Ultrasound therapy can be administered in two different ways: direct contact (mobile or fixed head technique) and immersion.

In the first mode, the area to be treated should be positioned so that it is relaxed, bare and dry. Apply therapeutic gel or cream.

The direct contact mode, which is most frequently used, consists in the application of the emitting head in direct contact with the skin with the interposition of a substance (usually a special conductive gel) to favor a on the side the transmission between the head and the skin and on the other the adherence, the slipping and the elimination of possible air between the skin and the transducer which could obstruct, due to its reflective capacity, the transmission of the ultrasonic wave. In the mobile head technique, the transducer, applied to the area to be treated with light pressure, is made to slide with circular movements. In the fixed head technique, the latter is placed on the area to be treated in contact with the skin for the entire duration of the therapy with the interposition of an electroconductive gel or therapeutic cream. With this technique, rapid increases in temperature are obtained in a very limited area, while the remaining field is not heated; therefore a lower supply power or a pulsed emission is required.

Immersion method

The part to be treated is immersed in a container containing water at a temperature of about 37 ° together with the emitting head, placed at a maximum distance of 2-3 cm from the body surface.

Caution

The parameters defined in the treatment programs (time, power and duty cycle) are intended to be applied only in the mobile head and immersion technique, where it is specified. In the fixed head technique, they must be considered reduced by 50%

Therapeutic indications

Not being dangerous, this therapy can also be applied to the elderly and children; the only contraindications concern PACE-MAKER carriers, heart failure patients, pregnant women. Below is a list of the main indications:

TRAUMATIC INJURIES It is used to relieve pain and to accelerate the healing of all fractures, sprains and dislocations caused by accidents, sports activities, including lumbar trauma.

ARTHRO-RHEUMATIC DISEASES It is used to treat myalgia, lumbago, arthrosis, epicondylitis and all bone pains of different nature, localized to the knee, hips, elbow, heel, hands, etc.

BONE CONSOLIDATION Accelerates the bone recalcification process in those who have suffered fractures and acts effectively in the treatment of femoral osteotomy, patella fractures, knee or arm surgery, neck or vertebrae injuries.

REGENERATION OF THE BONE From the analysis of radiographic plates of those who have undergone surgical interventions following traffic accidents or at work, it was found that with magnetotherapy, a process of regeneration and welding of the bone tissue was determined.

TENDINITIS It is used to treat all forms of ligament sprains of the knee, elbow, wrist, etc.

Contraindications

It is not recommended to use in the presence of metal synthesis means and joint prostheses, due to their greater absorption compared to the surrounding tissues and therefore for their possible deterioration or damage to neighboring tissues. In addition, localized treatment on the implant site can lead to overheating or to one

Contraindications

It is not recommended to use in the presence of metal synthesis means and joint prostheses, due to their greater absorption compared to the surrounding tissues and therefore for their possible deterioration or damage to neighboring tissues. In addition, localized treatment on the implant site can lead to overheating or to one

The absolute contraindication in the presence of osteoporosis appears to have fallen as the piezoelectric effect could even have a positive action on osteogenesis. It is contraindicated in case of acute inflammatory processes, neoplasms, skin lesions and sensitivity alterations. It should be avoided in correspondence with the eyeballs and the pregnancy uterus due to the possibility that the effect of cavitation may occur even at therapeutic dosages. Particular caution should be taken when treating areas near the spinal cord after laminectomy.

Therapy should be avoided in the cardiac area due to the possibility of interference on cardiac conduction and contraction; direct exposure of pacemakers or other implanted electromechanical equipment should be avoided due to the possible permanent damage that may result.

Specific fitness and beauty applications

- Applications for the absorption of active substances on the face and body The program is used for the application of cosmetic agents. Ultrasound facilitates the absorption of oils, fat-soluble vitamins, liposome products, emulsions and water-soluble agents.
- The effect is to increase cell permeability, decrease the functional barrier of the skin and increase the activity of the channels of the active ingredients.
- Treatment of superficial wrinkles and sensitive skin
- The treatment is used to stretch-smooth wrinkles [expression lines]. Thanks to its thermal effects, ultrasound heats the tissues and promotes biochemical and metabolic processes.
- Treatment of acne, pimples on the face and pimples on the body
- In addition to an anti-inflammatory effect, ultrasound activates the fibroblasts which play an important role in the post-inflammatory regeneration process. The collagen and elastic fibers produced after ultrasound treatment increase the elasticity of the tissues, reducing the possible formation of post acne scars. It is also possible to promote the absorption of anti-acne products by increasing their effect.
- **Anti cellulite treatments**
- Ideal for cellulite problems. There is a mobilization action of the fat in the treated area thanks to the thermal effect. The metabolism of fat and its removal is thus favored.

Specific health applications

The therapeutic effects of ultrasound are partly due to the increase in temperature. They are represented by analgesia with a fibrolytic and trophic effect.

- a) **Analgesia**
- The analgesic effect is due to the action of heat and probably also to a direct action of ultrasound on sensitive nerve endings.
- b) **Fibrolytic action**
- The oscillations of the particles of the tissues, produced by the ultrasounds, determine the disruption of the collagen fibers of the fibrous tissues.
- c) **Trophic effect**
- Vasodilation, which follows the thermal elevation, facilitates the removal of catabolites and delivers nutrients and oxygen to the tissues; in this way, ultrasounds improve tissue trophism, facilitate the repair of tissue damage and accelerate the resolution of inflammatory processes.


Sport specific applications

- **Release of contracted muscles**
- The relaxation of contracted muscles is linked to the thermal effect and the action of tissue micro-massage induced by ultrasound. Ideal for relaxing stressed muscles and eliminating muscle spasms caused by sports or injuries.

PATHOLOGY

TEMPOROMANDIBULAR PAIN


FREQUENCY	3 MHz
POWER	0.5 – 1 W/cm ²
EMISSION	puls. 50%
DURATION max	11 minutes
PROGRAM	53

PATHOLOGY

PATHOLOGIES DELTOID MUSCLE


FREQUENCY	1 MHz
POWER	0.5 – 1 W/cm ²
EMISSION	continuous
DURATION max	10 minutes
PROGRAM	59

PATHOLOGY

OCCIPITAL NEURALGIA


FREQUENCY	3 MHz
POWER	0.5 – 1 W/cm ²
EMISSION	puls. 75%
DURATION max	11 minutes
PROGRAM	54

PATHOLOGY

ROTATOR CUFF. MUSCLE DISTORSION


FREQUENCY	1 MHz
POWER	0.5 – 1 W/cm ²
EMISSION	puls. 80-100%
DURATION max	12 minutes
PROGRAM	60

PATHOLOGY

MYALGIA TRAPEZIUM


FREQUENCY	1 MHz
POWER	0.5 – 1 W/cm ²
EMISSION	puls. 50%
DURATION max	12 minutes
PROGRAM	55

PATHOLOGY

EPICONDYLITIS


FREQUENCY	3 MHz
POWER	0.5 – 1 W/cm ²
EMISSION	puls. 80-100%
DURATION max	10 minutes
PROGRAM	61

PATHOLOGY

PAIN AROUND SHOUL. BLADE


FREQUENCY	1 MHz
POWER	0.5 – 1 W/cm ²
EMISSION	puls. 75%
DURATION max	12 minutes
PROGRAM	56

PATHOLOGY

EPITROCLEITIS


FREQUENCY	3 MHz
POWER	0.5 – 1 W/cm ²
EMISSION	puls. 85%
DURATION max	10 minutes
PROGRAM	62

PATHOLOGY

TENDINITIS OF THE BICEPS

FREQUENCY	1 MHz
POWER	0.5 – 1 W/cm ²
EMISSION	continuous
DURATION max	12 minutes
PROGRAM	57

TREATMENT OF PATHOLOGIES WITH SONOPHORESIS

FREQUENCY	1 MHz
POWER	0.5 – 1 W/cm ²
EMISSION	continuous
DURATION max	10 minutes
PROGRAM	63


1 - APPLY THE SPECIFIC PRODUCT ON THE BODY IN REPLACEMENT OF ULTRASONIC GEL

2 - ACT WITH THE HANDPIECE ON THE PART OF THE BODY TO BE TREATED WITH CIRCULAR MOVEMENTS

PATHOLOGY

CAPSULITES-BORSITIS SHOUL. JOINTS


FREQUENCY	1 MHz
POWER	0.5 – 1 W/cm ²
EMISSION	puls. 50%
DURATION max	10 minutes
PROGRAM	58

PATHOLOGY

TENDINOPATHY BRACHIAL BICEPS

FREQUENCY	1 MHz
POWER	0.5 – 1 W/cm ²
EMISSION	puls. 90%
DURATION max	10 minutes
PROGRAM	64

PATHOLOGY

TENDINOPATHY OF THE BRACHIAL TRICEPS


FREQUENCY	3 MHz
POWER	0.5 – 1 W/cm ²
EMISSION	continuous
DURATION max	10 minutes
PROGRAM	65

If chronic

FREQUENCY	1 MHz
POWER	0.5 – 1 W/cm ²
EMISSION	pulsata al 50%
DURATION max	10 minutes
PROGRAM	66

Immersion (NO DIFFUSER)

FREQUENCY	1 MHz
POWER	2 – 2,5 W/cm ²
EMISSION	continuous
DURATION max	12 minutes
PROGRAM	67


PATHOLOGY

POLYARTHRITIS PRIMARY

Immersion (NO DIFFUSER)

FREQUENCY	1 MHz
POWER	1,5 - 2 W/cm ²
EMISSION	continuous
DURATION max	12 minutes
PROGRAM	72


PATHOLOGY

TRAPEZIOMETACARPAL METAC. ARTHROSIS

Immersion (NO DIFFUSER)


FREQUENCY	1 MHz
POWER	1,5 W/cm ²
EMISSION	continuous
DURATION max	12 minutes
PROGRAM	73

PATHOLOGY

MYOTENDONIC PATHOLOGIES FOREARM


FREQUENCY	3 MHz
POWER	0.5 – 1 W/cm ²
EMISSION	puls. 70-90%
DURATION max	10 minutes
PROGRAM	68

PATHOLOGY

DUPUYTREN'S DISEASE

FREQUENCY	3 MHz
POWER	0,8 - 1,2 W/cm ²
EMISSION	puls. 50%
DURATION max	11 minutes
PROGRAM	74


PATHOLOGY

TENDINITOPATHY HAND EXTENSORS

Immersion (NO DIFFUSER)


FREQUENCY	1 MHz
POWER	2 – 2,5 W/cm ²
EMISSION	puls. 85%
DURATION max	12 minutes
PROGRAM	69

PATHOLOGY

LOMBAGO

FREQUENCY	1 MHz
POWER	1 W/cm ²
EMISSION	puls. 70-90%
DURATION max	12 minutes
PROGRAM	75


PATHOLOGY

DE QUERVAIN TENDINOPATHY

Immersion (NO DIFFUSER)


FREQUENCY	1 MHz
POWER	2 W/cm ²
EMISSION	continuous
DURATION max	10 minutes
PROGRAM	70

PATHOLOGY

SACROILIAC PAIN

FREQUENCY	1 MHz
POWER	0.8 – 1 W/cm ²
EMISSION	puls. 80-100%
DURATION max	10 minutes
PROGRAM	76


PATHOLOGY

HAND FLEXOR TENDINITIS CARPAL TUNNEL SYNDROME

Immersion (NO DIFFUSER)


FREQUENCY	1 MHz
POWER	1,5-2 W/cm ²
EMISSION	puls. 85%
DURATION max	12 minutes
PROGRAM	71

PATHOLOGY

TENDINITIS ANTERIOR RECTUS


FREQUENCY	3 MHz
POWER	0.8 – 1 W/cm ²
EMISSION	puls. 40-60%
DURATION max	10 minutes
PROGRAM	77

PATHOLOGY

COXARTHROSIS PAIN


FREQUENCY	1 MHz
POWER	1 - 1,5 W/cm ²
EMISSION	puls. 75%
DURATION max	12 minutes
PROGRAM	78

PATHOLOGY

TENDINITIS - BUNS OF THE GOOSE LEG


FREQUENCY	3 MHz
POWER	0,8 - 1,2 W/cm ²
EMISSION	puls. 75%
DURATION max	10 minutes
PROGRAM	84

PATHOLOGY

PATHOLOGIES QUADRICEPS MUSCLE


FREQUENCY	1 MHz
POWER	0,8 - 1,5 W/cm ²
EMISSION	puls. 75%
DURATION max	12 minutes
PROGRAM	79

PATHOLOGY

FEMURAL-TIBIAL PATHOLOGIES


FREQUENCY	1 MHz
POWER	1 - 1,5 W/cm ²
EMISSION	puls. 75%
DURATION max	10 minutes
PROGRAM	85

PATHOLOGY

PATHOLOGIES FLEXOR MUSCLES


FREQUENCY	1 MHz
POWER	0,5 - 1,5 W/cm ²
EMISSION	puls. 75%
DURATION max	10 minutes
PROGRAM	80

PATHOLOGY

FEMURAL-TIBIAL LIGAMENT PATHOLOGIES


FREQUENCY	1 MHz
POWER	1 - 1,5 W/cm ²
EMISSION	puls. 50-70%
DURATION max	10 minutes
PROGRAM	86

PATHOLOGY

TENDINITIS PATELLA


FREQUENCY	3 MHz
POWER	0,8 - 1,2 W/cm ²
EMISSION	puls. 50-70%
DURATION max	10 minutes
PROGRAM	81

PATHOLOGY

FEMURAL BICEPS TENDINOPATHY


FREQUENCY	1 MHz
POWER	0,5 - 1 W/cm ²
EMISSION	puls. 75%
DURATION max	12 minutes
PROGRAM	87

PATHOLOGY

TENDINITIS QUADRICEPS


FREQUENCY	3 MHz
POWER	0,5 - 1,5 W/cm ²
EMISSION	puls. 50-60%
DURATION max	11 minutes
PROGRAM	82

PATHOLOGY

TENDINOPATHY POPLITHEUM


FREQUENCY	1 MHz
POWER	0,5 - 1 W/cm ²
EMISSION	puls. 50%
DURATION max	10 minutes
PROGRAM	88

PATHOLOGY

FEMURAL-PATELLA CHONDROPATHY


FREQUENCY	1 MHz
POWER	1 - 1,5 W/cm ²
EMISSION	continuous
DURATION max	10 minutes
PROGRAM	83

PATHOLOGY

SURAL TRICEPS PATHOLOGIES

FREQUENCY	1 MHz
POWER	0,8 - 1,2 W/cm ²
EMISSION	puls. 80%
DURATION max	12 minutes
PROGRAM	89


PATHOLOGY

TENDINITIS OF THE ANTERIOR TIBIAL

With sonophoresis

FREQUENCY	1 MHz
POWER	0,5 – 1,2 W/cm ²
EMISSION	continuous
DURATION max	10 minutes
PROGRAM	90


PATHOLOGY

ACHILLO-TIBIAL TENDINOPATHY

Immersion (NO DIFFUSER)


FREQUENCY	1 MHz
POWER	2 W/cm ²
EMISSION	puls. 75-100%
DURATION max	12 minutes
PROGRAM	94

PATHOLOGY

TENDINITIS OF FIBULA

FREQUENCY	1 MHz
POWER	0,5 – 1 W/cm ²
EMISSION	puls. 70%
DURATION max	12 minutes
PROGRAM	91


PATHOLOGY

HALUX VARUS TENDINOPATHY

Immersion (NO DIFFUSER)

FREQUENCY	1 MHz
POWER	2 W/cm ²
EMISSION	continuous
DURATION max	10 minutes
PROGRAM	95


PATHOLOGY

LATERAL ANKLE

By immersion (if edema or blood is present but in the acute and subacute phase) NO DIFFUSER

FREQUENCY	1 MHz
POWER	2 W/cm ²
EMISSION	puls. 50-70%
DURATION max	11 minutes
PROGRAM	92


PATHOLOGY

PLANTAR FASCIITIS

Immersion (NO DIFFUSER)

FREQUENCY	1 MHz
POWER	2 - 2,5 W/cm ²
EMISSION	puls. 50%
DURATION max	12 minutes
PROGRAM	96

PATHOLOGY

RETRO-ACHILLO BURSITIS

Immersion (NO DIFFUSER)

FREQUENCY	1 MHz
POWER	2 W/cm ²
EMISSION	continuous
DURATION max	10 minutes
PROGRAM	93


PATHOLOGY

METATARSALGIA METATARS.

Immersion (NO DIFFUSER)

FREQUENCY	1 MHz
POWER	2 - 2,5 W/cm ²
EMISSION	puls. 50-70%
DURATION max	10 minutes
PROGRAM	97

PROBE TESTING

2-minute program to be used only and exclusively to check the correct functioning of the treatment handpiece

PROGRAM	98
---------	----

Questions and answers

Information

When can the appliance be used?

The device can be used at any time of the day, not only as a treatment against cellulite, but also as a prevention or maintenance of the results obtained. We recommend that you carve out a moment of the day, to be able to devote yourself exclusively to the treatment, without having to deal with anything else.

How long does it take to get results?

It depends; in most cases, even from the first applications. However, it must be considered that the individual's reactivity is subjective and varies from person to person, from the complexity of the imperfection and from the affected surface. It is advisable not to desist from applications and to provide for a mild and continuous physical activity and adequate nutrition, excluding drastic diets.

Can ultrasound be harmful if used often?

Caution

Ultrasound applications do not create any damage except in the contraindicated cases (see the chapter "Contraindications"). However, it is advisable to adhere to what is indicated, avoiding excessive actions or too frequent applications that could be aesthetically counterproductive.

If you increase the potency, do you get more results?

Caution

NO. The appliance is equipped with specific programs that already provide for the optimal power level to be used, the presence of buttons to decrease or increase allows for personal adjustment, but in no case does it allow to exceed the predefined power levels. Personal adjustment must include a pleasant and comfortable feeling, like a massage.

How long does it take to see the first results?

The appliance is equipped with an exclusive function, which allows you to act on any type of problem, with maximum effectiveness. The results depend on the stage of advancement of the blemish, however by selecting the most suitable treatment with the keyboard, you will be able to appreciate a first result after application, a smoother and more velvety skin to the touch.

The subsequent results, adopting the appropriate precautions described in the manual, will soon be seen, bearing in mind however that cellulite, in some cases, requires a minimum of two complete cycles of applications per year (spring and autumn). . It is advisable to drink plenty of water, possibly non-carbonated, to adopt an appropriate diet and a mild but continuous physical activity (walking and cycling) avoiding both drastic diets and excessive and too prolonged physical exertion.

If you do not feel anything during the treatment or you do not perceive a continuous and light vibration, what does this mean?

It means that the area in question is less sensitive, but this does not mean that the selective ultrasound wave of the device is not performing its function. The selectivity of the device's system acts at different depths and varies continuously, so your sensation will also vary. Trust the chosen program and you will get the best possible result.

When can the maintenance program be used?

It is advisable to use the maintenance program alternating it with every three applications of specific treatment, until satisfactory results are obtained. Later you can use this program - but whenever you want.

Is it possible to do more treatments on the same area several times a day?

What is the maximum recommended time of use per day?

We recommend daily use of a maximum duration of 15 minutes for each body area (programs with a duration of 6 minutes can also be repeated twice).

Can I hold the treatment head still?

Caution

NO, if the head is not slid, the part being treated could overheat, both annoying and having a negative effect. Therefore it is necessary to follow the method of application indicated, to carry out safe treatments, otherwise the effectiveness of the ultrasound action would be limited to a single point, while it is advisable to act on the whole surrounding area, reactivating the circulation on the whole the treated area.

Should I always apply the gel or cream to the treatment head and skin?

Caution

Yes, and in sufficient quantity for the entire area affected by the application, otherwise the ultrasounds are not transmitted and the treatment head could heat up too much and not exert its aesthetic effect and more importantly create discomfort on the skin.

How do I know if the head is working properly?

If you have doubts about the correct operation and production of ultrasounds (sounds that are not audible to the human ear) select program 98 called "**V-SONO PLUS**" and let a few drops of water fall on the head which will appear in motion, it will seem to fry as well as presenting a slight evaporation if the large handpiece is inserted. If nothing happens, see the "**Troubleshooting**" chapter.

How can I clean the treatment head after each application?

Cleaning the treatment head is extremely important and must be carefully performed after each application. The residues of gels or creams made to dry can induce the birth of bacteria and molds, in addition to limiting the correct functioning of the ultrasonic transducer incorporated in the head. For cleaning, follow the instructions in the "Cleaning" chapter.

If the treatment head has cracks or the cable appears damaged, can I still use them?

Caution

Absolutely no! After each treatment, the treatment head must be checked as indicated in the "**Maintenance**" chapter.

Can I have applications on the lower abdomen?

Caution

NO, these applications are absolutely contraindicated and not recommended; in this area you are close to noble organs and especially in the case of severe dance or menstruation you must absolutely not make applications of any kind (see the chapter "Contraindications").

Can applications be made on children or adolescents?

Caution

NOT, it is absolutely contraindicated to act with any means in the presence of growing cartilage (see the "**Contraindications**" chapter).

Can I do applications on varicose veins?

Caution

NO, any type of treatment is absolutely contraindicated if you have severe circulatory problems and in the presence of such pathologies. However, it can treat other areas of the body not directly affected. In the most serious cases it is advisable to consult a doctor (see the "**Contraindications**" chapter).

I have a pacemaker, can I use the device?

Caution

NO the use of electronic devices in general is contraindicated for people with cardiac pacemakers and also in the presence of implanted electrodes. It is advisable to consult a doctor (see the "**Contraindications**" chapter).

Troubleshooting

⚠ Caution

No intervention on the equipment is permitted, as in addition to generating possible dangerous situations, it will invalidate the validity of the Warranty.

What to do if the appliance does not show any sign of operation when the power switch is turned on?

Check that the power cable is properly inserted in its housing and in the socket.

Check the power cable for signs such as:

- - breaks or interruptions;
- deterioration.

Check the integrity of the fuse inside the fuse holder.

⚠ Caution

In case of breakage or deterioration of the power cable, immediately disconnect it from the power socket and replace it immediately.

After having turned on the general switch the appliance does not enable the functions?

Check the treatment head connection.

Screw in the treatment head connector.

Select programme N. 98 (**V-SONO PLUS**).

Turn on the **START** key and if no effect is noted, check that the **"PAUSE"** and **"STOP"** led are off and that the **"ACTIVE"** led is on.

If the **"ACTIVE"** led is off, repeat the procedure.

If the **"ACTIVE"** led is on, check the emission of ultrasounds by allowing a few drops of water to fall on the head; it will appear to move, seem to fry and even evaporate slightly if the large knob is plugged in (see photo below).

Reference regulations

Appliance has been produced in compliance with the following regulations:

- Directive 93/42/EEC, dated 14th June 1993, regarding medical appliances, adopted by Legislative Decree no. 46 dated 24/02/97 and its subsequent modifications.

CEI EN 60601-1
CEI EN 60601-1-2
CEI EN 60601-2-5
CEI EN 60601-1-11
CEI EN 14971

Declaration of conformity

The company HSD SRL, with registered and operational headquarters in SERRAVALLE (RSM) in Strada Acquisalata, n 11 internal 6, declares, under its sole responsibility, that the medical device called risk class IIa in accordance with rule 9 of Annex IX of Directive 93/42 / EEC and subsequent amendments (implemented in Italy with Legislative Decree No. 46 of 24 February 1997, and subsequent amendments)

- complies with the essential requirements and provisions of Directive 93/42 / EEC and subsequent amendments and additions. as per the Technical File filed at the Company;

- is manufactured in accordance with the Quality System that meets the requirements set out in Annex V + VII of the aforementioned Legislative Decree, as per the Certificate issued by BUREAU VERITAS ITALIA S.p.a. - Notified Body n. 1370
- Viale Monza, n. 347 - 20126 Milan .

- complies with Directive 2011/65 / EU of the European Parliament and of the Council of 8 June 2011, on the restriction of the use of certain hazardous substances in electrical and electronic equipment.

- the Authorized European Representative is Ms Manuela Corbelli Via Boccioni n. 20- Morciano di Romagna

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Serravalle (RSM)

NOTE

Lined area for notes, consisting of 30 horizontal lines.

NOTE

Lined area for notes, consisting of 30 horizontal lines.

Fabbricante: HSD S.r.l.
Strada Acquasalata, 11 - int. 6 - 47899
Serravalle RSM
T. (+378) 0549.941496 - F. (+378) 0549.976456

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CAUTION: please keep the purchase document (receipt or invoice) in order to confirm the purchase date and exercise the guarantee rights, in compliance with existing legislation.

Guarantee certificate *Part for the customer to fill in and send in a sealed envelope*

Code

Model

Serial no.

Surname

Name

Street

No.

City

Postcode

Country

Telephone no.

Fault found:

Guarantee certificate *Part for the customer to fill in and keep*

Code

Model

Serial no.

Surname

Name

Street

No.

City

Postcode

Country

Telephone no.

MADE IN ITALY

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© 2019 HSD S.r.l. RSM - Edizione Luglio 2019 - Rel. 0.1